



# SAVVY<sup>®</sup> & SAVVY<sup>®</sup> Long



Versatile 0.018"  
over-the-wire balloons  
designed to meet  
the challenges of  
femoropopliteal and  
infrapopliteal lesions



# SAVVY<sup>®</sup> & SAVVY<sup>®</sup> Long

## Versatile 0.018" over-the-wire balloons designed to meet the challenges of femoropopliteal and infrapopliteal lesions

### Ordering Information

#### SAVVY<sup>®</sup>

Balloon OD x Length (mm x cm)	Rated Burst Pressure (atm)	Recommended Sheath Fit (F)	Ordering Number	Shaft Length (cm)		
				80	120	150
2 x 2	10	4	435-202	S	L	-
2 x 4	10	4	435-204	S	L	X
2 x 6	10	5	435-206	S	L	-
2 x 10	10	5	435-200	S	L	X
2.5 x 2	10	4	435-252	S	L	-
2.5 x 4	10	4	435-254	S	L	X
2.5 x 6	10	5	435-256	S	L	-
2.5 x 10	10	5	435-250	S	L	-
3 x 2	10	4	435-302	S	L	-
3 x 4	10	4	435-304	S	L	X
3 x 6	10	5	435-306	S	L	-
3 x 10	10	5	435-300	S	L	X
3.5 x 2	10	4	435-352	S	L	-
3.5 x 4	10	4	435-354	S	L	X
3.5 x 6	10	5	435-356	S	L	-
3.5 x 10	10	5	435-350	S	L	-
4 x 2	10	4	435-402	S	L	-
4 x 4	10	4	435-404	S	L	X
4 x 6	10	5	435-406	S	L	-
4 x 10	10	5	435-400	S	L	X
4.5 x 2	10	4	435-452	S	L	-
4.5 x 4	10	4	435-454	S	L	-
4.5 x 6	10	5	435-456	S	L	-
5 x 2	10	4	435-502	S	L	-
5 x 4	10	4	435-504	S	L	-
5 x 6	10	5	435-506	S	L	-
5.5 x 2	10	4	435-552	S	L	-
5.5 x 4	10	4	435-554	S	L	-
5.5 x 6	10	5	435-556		L	-
5 x 10	10	5	435-500	S	L	-
6 x 2	10	4	435-602	S	L	-
6 x 4	10	4	435-604	S	L	-
6 x 6	10	5	435-606	S	L	-
6 x 10	10	5	435-600	S	L	-

#### SAVVY<sup>®</sup> Long

Balloon OD x Length (mm x cm)	Rated Burst Pressure (atm)	Recommended Sheath Fit (F)	Ordering Number	Shaft Length (cm)		
				80	120	150
2 x 12	15	4	436-2012	S	L	-
2 x 15	15	4	436-2015	S	L	-
2 x 22	15	4	436-2022	S	L	-
2.5 x 12	15	4	436-2512	S	L	-
2.5 x 15	15	4	436-2515	S	L	-
2.5 x 22	15	4	436-2522	S	L	-
3 x 12	15	4	436-3012	S	L	-
3 x 15	15	4	436-3015	S	L	-
3 x 22	15	4	436-3022	S	L	-
3.5 x 12	15	4	436-3512	S	L	-
3.5 x 15	15	4	436-3515	S	L	-
3.5 x 22	15	4	436-3522	S	L	-
4 x 12	15	4	436-4012	S	L	-
4 x 15	15	4	436-4015	S	L	-
4 x 22	15	4	436-4022	S	L	-
5 x 12	13	4	436-5012	S	L	-
5 x 15	13	4	436-5015	S	L	-
5 x 22	13	4	436-5022	S	L	-
6 x 12	12	5	436-6012	S	L	-
6 x 15	12	5	436-6015	S	L	-
6 x 22	12	5	436-6022	S	L	-

