



# SLEEK<sup>™</sup>

The first and only 0.014"  
rapid-exchange balloon  
fully dedicated to  
infrapopliteal lesions



# SLEEK™

The first and only 0.014" rapid-exchange balloon fully dedicated to infrapopliteal lesions

## Ordering Information

Balloon OD x Length (mm x cm)	Rate Burst Pressure (atm)	Recommended Sheath Fit (F)	Ordering Number	Shaft Length 150 (cm)
2 X 4	16	4	425-2004	X
2 x 6	16	4	425-2006	X
2 x 12	15	4	425-2012	X
2 x 15	15	4	425-2015	X
2 x 22	15	4	425-2022	X
2.5 x 4	16	4	425-2504	X
2.5 x 6	16	4	425-2506	X
2.5 x 12	15	4	425-2512	X
2.5 x 15	15	4	425-2515	X
2.5 x 22	15	4	425-2522	X
3 x 4	16	4	425-3004	X
3 x 6	16	4	425-3006	X
3 x 12	15	4	425-3012	X
3 x 15	15	4	425-3015	X
3 x 22	15	4	425-3022	X
3.5 x 4	16	4	425-3504	X
3.5 x 6	16	4	425-3506	X
3.5 x 12	15	4	425-3512	X
3.5 x 15	15	4	425-3515	X
3.5 x 22	15	4	425-3522	X
4 x 4	16	4	425-4004	X
4 x 6	16	4	425-4006	X
4 x 12	15	4	425-4012	X
4 x 15	15	4	425-4015	X
4 x 22	15	4	425-4022	X

