

UNMOUNTED PALMAZ GENESIS STENT (MEDIUM)

ORDERING INFORMATION

UNMOUNTED PALMAZ GENESIS (MEDIUM)

EXPANSION RANGE FROM 5 TO 8 MM

<i>Unexpanded Stent Length (mm)</i>	<i>Expanded Stent Length (mm)</i>	<i>Expanded Stent diameter (mm)</i>	<i>Guiding Compatibility 8F / 1.088"</i>	<i>Sheath Compatibility (F)</i>	<i>Catalogue Number</i>
12	12	5	8	6	PG124P
12	11	6	8	6	
12	10	7	8	7	
12	10	8	8	7	
15	15	5	8	6	PG154P
15	14	6	8	6	
15	13	7	8	7	
15	13	8	8	7	
18	18	5	8	6	PG184P
18	17	6	8	6	
18	17	7	8	7	
18	15	8	8	7	
24	24	5	8	6	PG244P
24	23	6	8	6	
24	23	7	8	7	
24	21	8	8	7	

Recommended for mounting on Cordis Opta™ Pro PTA Balloon Catheter.

COMPLEMENTARY CORDIS PRODUCTS

In order to further enhance its procedural efficiency and performance, we recommend to combine Palmaz® Genesis™ with other Cordis product families. The products listed hereunder are only a selection. Please contact your local sales representative for further advice on matching products.

ORDERING INFORMATION

JINDO™


STEERABLE GUIDEWIRES



<i>Length (cm)</i>	<i>Taper</i>	<i>Catalogue Numbers</i>
300	Short	503-453
300	Intermediate	503-553

VISTA BRITE TIP® IG

INTRODUCER GUIDES



<i>Catheter Shape</i>	<i>Length (cm)</i>	<i>OD (F)</i>	<i>ID</i>	<i>Catalogue Number</i>
Renal curve (small)	55	8	.088"	403-8553R
Renal curve (adult)	55	8	.088"	403-8553A

Important Information: Prior to use, refer to the "Instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

Europe
Tel: +32 2 352 14 11
Fax: +32 2 352 15 90