

UNMOUNTED PALMAZ GENESIS STENT (LARGE)

ORDERING INFORMATION

UNMOUNTED PALMAZ GENESIS (LARGE)

EXPANSION RANGE FROM 5 TO 10 MM

Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Expanded stent diameter (mm)	Guiding Compatibility (F)	Sheath Compatibility (F)	Catalogue Number
29	28	5	8 (.088")	7	PG295P
29	28	6	9 (.098")	7	
29	27	7	9 (.098")	7	
29	26	8	9 (.098")	7	
29	25	9	10 (.110")	7.5	
29	24	10	10 (.110")	7.5	
39	38	5	8 (.088")	7	PG395P
39	38	6	9 (.098")	7	
39	37	7	9 (.098")	7	
39	36	8	9 (.098")	7	
39	35	9	10 (.110")	7.5	
39	32	10	10 (.110")	7.5	
59	59	5	8 (.088")	7	PG595P
59	58	6	9 (.098")	7	
59	57	7	9 (.098")	7	
59	56	8	9 (.098")	7	
59	53	9	10 (.110")	7.5	
59	53	10	10 (.110")	7.5	
79	79	5	8 (.088")	7	PG795P
79	78	6	9 (.098")	7	
79	76	7	9 (.098")	7	
79	76	8	9 (.098")	7	
79	75	9	10 (.110")	7.5	
79	73	10	10 (.110")	7.5	

EXPANSION RANGE FROM 10 TO 12 MM

Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Expanded stent diameter (mm)	Guiding Compatibility (F)	Sheath Compatibility (F)	Catalogue Number
19	17	10	10 (.110")	8	PG1910P
19	16	12	10 (.110")	8	
29	27	10	10 (.110")	8	PG2910P
29	25	12	10 (.110")	8	
39	36	10	10 (.110")	8	PG3910P
39	34	12	10 (.110")	8	
59	55	10	10 (.110")	8	PG5910P

Recommended for mounting on Cordis Opta™ Pro PTA Balloon Catheter.

COMPLEMENTARY CORDIS PRODUCTS

In order to further enhance its procedural efficiency and performance, we recommend to combine Palmaz® Genesis™ with other Cordis product families. The products listed hereunder are only a selection. Please contact your local sales representative for further advice on matching products.


ORDERING INFORMATION

.035" STORQ™ STEERABLE GUIDEWIRES



<i>Diameter (in)</i>	<i>Tip Flexibility</i>	<i>Tip Shape</i>	<i>Length (cm)</i>	<i>Catalogue Number</i>
.035	Standard	Straight	180	503-456
.035	Standard	Angled	180	503-456J
.035	Standard	Modified J	180	503-456MJ
.035	Standard	Straight	300	503-456X
.035	Standard	Angled	300	503-456Y
.035	Standard	Modified J	300	503-456MY

BRITE TIP® SHEATH INTRODUCERS



<i>Cannula usable length (cm)</i>	<i>Catalogue Number 6F</i>	<i>Catalogue Number 7F</i>	<i>Catalogue Number 8F</i>
11	401-611 M	401-711 M	401-811 M
23	401-623 M	401-723 M	401-823 M
35	401-635 M	401-735 M	401-835 M
45	401-645 M	401-745 M	401-845 M
55	401-655 M	401-755 M	401-855 M
90	401-690 M	401-790 M	401-890 M

Important Information: Prior to use, refer to the "Instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

Europe
Tel: +32 2 352 14 11
Fax: +32 2 352 15 90