


Cordis Long S.M.A.R.T.[®]

new
120 & 150
mm
LENGTH

Think S.M.A.R.T.[®]



ALL LENGTHS COMPATIBLE WITH 6F SHEATH



Introducing the longer S.M.A.R.T.®

Since their introduction in 1999, S.M.A.R.T.® Stents have built and expanded their leadership and their legacy. Every S.M.A.R.T.® Stent is a reaffirmation of the Cordis Endovascular commitment to better patient care through evidence-based, results-oriented innovation.



WE'RE NOT JUST INCREASING THE LENGTH.


WE'RE GROWING THE LEGACY.



A legacy of proven design¹

- **SUPERB SCAFFOLDING**
for uniform coverage
- **OFFSET ALTERNATING BRIDGES**
for wall apposition on a bend
- **SHORT STRUT LENGTH**
for individual cell response

¹ Cordis data on file



A legacy of confidence

- **MORE THAN 800,000 STENT PLACEMENTS WORLDWIDE**
- **OUTSTANDING CLINICAL EVIDENCE⁽¹⁻²⁾**
- **THE MOST EXTENSIVELY STUDIED NITINOL STENT⁽¹⁻²⁾**

¹SIROCCO II, J. Endovasc. Ther. 2006; 13:701-710
²CRISP_US, J. Vasc. Interv. Radiol. 2004; 15:911-918

Think **S.M.A.R.T.®**

Ordering Information

Cordis S.M.A.R.T.® Stents Long 120mm and 150mm Stents



| 120cm Delivery System | Expanded Stent Diameter (mm) | Expanded Stent Length (mm) | Profile: Sheath/Guide Fit (F) |
|-----------------------|------------------------------|----------------------------|-------------------------------|
| C06120MV | 6 | 120 | 6/8 |
| C07120MV | 7 | 120 | 6/8 |
| C08120MV | 8 | 120 | 6/8 |
| C06150MV | 6 | 150 | 6/8 |
| C07150MV | 7 | 150 | 6/8 |
| C08150MV | 8 | 150 | 6/8 |

Cordis S.M.A.R.T.® CONTROL™ Stents: Large, 12 x 14mm diameter

| 80cm Delivery System | 120cm Delivery System | Expanded Stent Diameter (mm) | Expanded Stent Length (mm) | Profile: Sheath/Guide Fit (F) |
|----------------------|-----------------------|------------------------------|----------------------------|-------------------------------|
| C12030SV | C12030MV | 12 | 30 | 7/9 |
| C12040SV | C12040MV | 12 | 40 | 7/9 |
| C12060SV | C12060MV | 12 | 60 | 7/9 |
| C12080SV | C12080MV | 12 | 80 | 7/9 |
| C14030SV | C14030MV | 14 | 30 | 7/9 |
| C14040SV | C14040MV | 14 | 40 | 7/9 |
| C14060SV | C14060MV | 14 | 60 | 7/9 |
| C14080SV | C14080MV | 14 | 80 | 7/9 |

Cordis S.M.A.R.T.® CONTROL™

| 80cm Delivery System | 120cm Delivery System | Expanded Stent Diameter (mm) | Expanded Stent Length (mm) | Profile: Sheath/Guide Fit (F) |
|----------------------|-----------------------|------------------------------|----------------------------|-------------------------------|
| C06020SV | C06020MV | 6 | 20 | 6/8 |
| C06030SV | C06030MV | 6 | 30 | 6/8 |
| C06040SV | C06040MV | 6 | 40 | 6/8 |
| C06060SV | C06060MV | 6 | 60 | 6/8 |
| C06080SV | C06080MV | 6 | 80 | 6/8 |
| C06100SV | C06100MV | 6 | 100 | 6/8 |
| C07020SV | C07020MV | 7 | 20 | 6/8 |
| C07030SV | C07030MV | 7 | 30 | 6/8 |
| C07040SV | C07040MV | 7 | 40 | 6/8 |
| C07060SV | C07060MV | 7 | 60 | 6/8 |
| C07080SV | C07080MV | 7 | 80 | 6/8 |
| C07100SV | C07100MV | 7 | 100 | 6/8 |
| C08020SV | C08020MV | 8 | 20 | 6/8 |
| C08030SV | C08030MV | 8 | 30 | 6/8 |
| C08040SV | C08040MV | 8 | 40 | 6/8 |
| C08060SV | C08060MV | 8 | 60 | 6/8 |
| C08080SV | C08080MV | 8 | 80 | 6/8 |
| C08100SV | C08100MV | 8 | 100 | 6/8 |
| C09020SV | C09020MV | 9 | 20 | 6/8 |
| C09030SV | C09030MV | 9 | 30 | 6/8 |
| C09040SV | C09040MV | 9 | 40 | 6/8 |
| C09060SV | C09060MV | 9 | 60 | 6/8 |
| C09080SV | C09080MV | 9 | 80 | 6/8 |
| C10020SV | C10020MV | 10 | 20 | 6/8 |
| C10030SV | C10030MV | 10 | 30 | 6/8 |
| C10040SV | C10040MV | 10 | 40 | 6/8 |
| C10060SV | C10060MV | 10 | 60 | 6/8 |
| C10080SV | C10080MV | 10 | 80 | 6/8 |

Important Information:
Prior to use, refer to the "Instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

© Cordis Europe, a division of Johnson & Johnson Medical NV/SA • Tel: +32 2 352 14 11 • January 2007 • 2E-800-1584-1