


SUPER TORQUE MB™

5F MARKER BAND FLUSH CATHETERS

The Gauge
for
Success




Cordis
a Johnson & Johnson company
ENDOVASCULAR

FOR A LESS INVASIVE FUTURE

Matching the exacting standards of Cordis quality, the Super Torque MB™ markerband catheter adds the possibility of performing precise length measurements to the widely renowned angiographic qualities of the Super Torque diagnostic catheters.

The excellent visibility of the gold-alloy marker bands provides easy determination of the exact length of a (covered) stent for stenotic arteries or an AAA-device for aortic aneurisms.

The braided body provides you with a highly manoeuvrable device that is demanded by this type of challenging procedure.



ORDERING INFORMATION

Super Torque MB™

# of marker bands	Length (cm)	Max pressure (psi)	Guidewire compatibility	Catalogue number
10	110	1200	.038"	532-598A
20	110	1200	.038"	532-598B

Important Information: Prior to use, refer to the "Instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of the Cordis policy of continuous product development we reserve the right to change product specifications without prior notification.

Europe
Tel: +32 2 352 14 11
Fax: +32 2 352 15 90