


Aviator™ Plus



PTA Dilatation Catheter
– for speed and
efficiency in carotid
and renal peripheral
endovascular
interventions



Aviator™ Plus

PTA Dilatation Catheter – for speed and efficiency in carotid and renal peripheral endovascular interventions

Ordering Information

Balloon OD x Length (mm x cm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Fit (F)	Guiding Catheter Fit (F)	Ordering Number	Shaft Length (142 cm)
4 x 1.5	10	14	4	6	424-4015	W
4 x 2	10	14	4	6	424-4020	W
4 x 3	10	14	4	6	424-4030	W
4 x 4	10	14	4	6	424-4040	W
4.5 x 1.5	10	14	4	6	424-4515	W
4.5 x 2	10	14	4	6	424-4520	W
4.5 x 3	10	14	4	6	424-4530	W
4.5 x 4	10	14	4	6	424-4540	W
5 x 1.5	10	14	4	6	424-5015	W
5 x 2	10	14	4	6	424-5020	W
5 x 3	10	14	4	6	424-5030	W
5 x 4	10	14	4	6	424-5040	W
5.5 x 1.5	10	14	4	6	424-5515	W
5.5 x 2	10	14	4	6	424-5520	W
5.5 x 3	10	14	4	6	424-5530	W
5.5 x 4	10	14	4	6	424-5540	W
6 x 1.5	10	14	4	6	424-6015	W
6 x 2	10	14	4	6	424-6020	W
6 x 3	10	14	4	6	424-6030	W
6 x 4	10	14	4	6	424-6040	W
7 x 1.5	10	12	5	6	424-7015	W
7 x 2	10	12	5	6	424-7020	W
7 x 3	10	12	5	6	424-7030	W
7 x 4	10	12	5	6	424-7040	W

