

ENVOY®

GUIDING CATHETER




Today. Tomorrow. Envoy.


**Cordis**  
a Johnson & Johnson company  
NEUROVASCULAR

# ENVOY® Guiding Catheters

- **ENVOY** and **ENVOY XB** Guiding Catheters provide treatment options for procedural versatility
- Soft radiopaque tip enables less traumatic vessel placement
- Lubricious PTFE inner liner provides smooth device passage
- Nylon/polyurethane and stainless steel braid shaft materials




**MULTI-SEGMENT™ Construction** balances distal flexibility with reliable torque and support.



Newly designed Hybrid Braid delivers maximum flow rates while maintaining catheter support and stability.


## 6F Guiding Catheter Flow Rates



Enhanced Large Lumen (0.070 in.) of the **6F ENVOY®** enables exceptional flow rates around in-dwelling devices.

Data on file, Cordis Neurovascular, Inc.

## Shapes



## Ordering Information

Catalog Number	Usable Length (cm)	Shape
<b>ENVOY® 5F - 0.056 in. ID (1.4mm)</b>		
556-26090	90	Straight
556-26000	100	Straight
556-25690	90	Multipurpose C
556-25600	100	Multipurpose C
556-25890	90	Multipurpose D
556-25800	100	Multipurpose D
556-25090	90	Modified Cerebral (Burke)
556-25000	100	Modified Cerebral (Burke)
556-25490	90	Headhunter 1
556-25400	100	Headhunter 1
<b>ENVOY® 6F - 0.070 in. ID (1.8mm)</b>		
670-26090	90	Straight
670-26000	100	Straight
670-25690	90	Multipurpose C
670-25600	100	Multipurpose C
670-25890	90	Multipurpose D
670-25800	100	Multipurpose D
670-25090	90	Modified Cerebral (Burke)
670-25000	100	Modified Cerebral (Burke)
670-25490	90	Headhunter 1
670-25400	100	Headhunter 1
670-25290	90	Simmons 2
670-25200	100	Simmons 2

## ENVOY® XB 6F - 0.070 in. ID (1.8mm)

670-26090B	90	Straight
670-26000B	100	Straight
670-25690B	90	Multipurpose C
670-25600B	100	Multipurpose C
670-25890B	90	Multipurpose D
670-25800B	100	Multipurpose D
670-25090B	90	Modified Cerebral (Burke)
670-25000B	100	Modified Cerebral (Burke)
670-25490B	90	Headhunter 1
670-25400B	100	Headhunter 1
670-25290B	90	Simmons 2
670-25200B	100	Simmons 2

USA  
**Cordis Neurovascular, Inc.**  
 14000 NW 57th Court  
 Miami Lakes, FL 33014  
 For Customer Service Call 1-800-327-7714

Europe  
**Johnson & Johnson Medical NV/SA**  
 Tel: +32 (0)2 352 14 11  
 Fax: +32 (0)2 352 15 90

**Cordis**  
 a Johnson & Johnson company  
 NEUROVASCULAR